


Product Name: Recombinant Human STK31 Protein
Catalog #: PEH8992



Summary

Name	STK31_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXU1
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NDAGNLITFP KESLAVGDFN LGSNVSLEKI KQDQKLIEN EKLKTEKDAL LESYKALELK VEQIAQELQQ EKAAAVDLTN HLEYTLKTYI DTRMKNLAAK MEILKEMRHV DISVRFGKDL SDAIQVLDEG CFTTPASLNG LEIWAEYSL AQENIKTCEY VSEGNILIAQ RNEMQQKLYM SVEDFILEVD ESSLNKRLKT LQDLSVSLEA VYGQAKEGAN SDEILKKFYD WKCDKREEFT SVRSETDASL HRLVAWFQRT LKVFDLSVEG SLISEDAMDN IDEILEKTES SVCKELEIAL VDQGDADKEI ISNTYSQVLQ KIHSEERLIA TVQAKYKDSI EFKKQLIEYL NKSPSVDHLL SIKKTLKSLK ALLRWKLVEK SNLEESDDPD GSQIEKIKEE

image



Background

Alternative Names Serine/threonine-protein kinase 31 (EC 2.7.11.1) (Serine/threonine-protein

Product Name: Recombinant Human STK31 Protein
Catalog #: PEH8992



kinase NYD-SPK) (Sugen kinase 396) (SgK396)

Note

For Research Use Only , Not for Diagnostic Use