


Product Name: Recombinant Human SH2D4A Protein
Catalog #: PEH8995



Summary

Name	SH24A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H788
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MLKQILSEMY IDPDLLAELS EEQKQILFFK MREEQIRRWK EREAAMERKE SLPVKPRPKK ENGKSVHWKL GADKEVWVWV MGEHHLDPKY DVLCEIIAE RARLKAEQEA EEPRKTHSEE FTNSLKTKSQ YHDLQAPDNQ QTKDIWKKVA EKEELEQGSR PAPTLEEEKI RSLSSSRNI QQMLADSINR MKAYAFHQKK ESMKKKQDEE INQIEEERTK QICKSWKEDS EWQASLRKSK AADEKRRSLA KQAREDYKRL SLGAQKGRGG ERLQSPLRVP QKPERPPLPP KPQFLNSGAY PQKPLRNQGV VRTLSSSAQE DIIRWFKEEQ LPLRAGYQKT SDTIAPWFHG ILTLKKANEL LLSTGMPGSF LIRVSEIRIG YALSYLSEDG CKHFLIDASA DAYSFLGVDQ LQHATLADLV EYHKEEPITS LGKELLYPC GQQDQLPDYL ELFE

image



Background

Product Name: Recombinant Human SH2D4A Protein
Catalog #: PEH8995



Alternative Names

SH2 domain-containing protein 4A (Protein SH(2)A) (Protein phosphatase 1 regulatory subunit 38)

Note

For Research Use Only , Not for Diagnostic Use