


Product Name: Recombinant Human SHC3 Protein
Catalog #: PEH8997



Summary

| | |
|------------------------------|--|
| Name | SHC3_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q92529 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 28 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | PPGGFLDTRL KPRPHAPDTA QFAGKEQTY Y QGRHLGDTFG EDWQQTPLRQ GSSDIYSTPE GKLHVAPTGE APTYVNTQQI PPQAWPAAVS SAESSPRKDL FDMKPFEDAL KNQPLGPVLS KAASVECISP VSPRAPDAKM LEELQAETWY QGEMSRKEAE GLEKDGDFL VRKSTTNP GS FVLTGMHNGQ AKHLLLDVPE GTIRTKDRVF DSISHLINHH LESSLPIVSA GSELCLQQPV |

image



Background

| | |
|--------------------------|---|
| Alternative Names | SHC-transforming protein 3 (Neuronal Shc) (N-Shc) (Protein Rai) (SHC-transforming protein C) (Src homology 2 domain-containing-transforming protein C3) (SH2 domain protein C3) |
|--------------------------|---|

Product Name: Recombinant Human SHC3 Protein
Catalog #: PEH8997



Note

For Research Use Only , Not for Diagnostic Use