


Summary

Name	SRG2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	PODMP2
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEIIAEYDTQ VKEIRAQLTE QMKCLDQQCE LRVQLLDLQ DFFRKKAEIE MDYSRNLEKL AERFLAKTCS TKDQQFKKDQ NVLSPVNCWN LLLNQVKRES RDHTTSLDIY LNNIIPRFVQ VSEDSGRLFK KSKEVGQQLQ DDLMKVLNEL YSVMKTYHMY NADSISAQSK LKEAEKQEEK QIGKSVKQED RQTPRSPDST ANVRIEEKHV RRSSVKKIEK MKEKRQAKYT ENKLKAIKAR NEYLLALEAT NASVFYYIH DLSDLIDCCD LGYHASLNRA LRTFLSAELN LEQSKHEGLD AIENAVENTL ATSDKQRLME MYNNVFCPPM KFEFQPHMGD

image



Background

Alternative Names	SLIT-ROBO Rho GTPase-activating protein 2B (SLIT-ROBO Rho GTPase activating protein 2 pseudogene 2)
--------------------------	---

Product Name: Recombinant Human SRGAP2B Protein
Catalog #: PEH8998



Note

For Research Use Only , Not for Diagnostic Use