


Product Name: Recombinant Human STRBP Protein
Catalog #: PEH9005



Summary

Name	STRBP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96SI9
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LD SKAIDL MN ALMRLNQIRP GLQYKLLSQS GPVHAPVFTM SVDVDGTTYE ASGPSKKTAK LHVAVKVLQA MGYPTGFDAD IECMSSDEKS DNESKNETVS SNSSNNTGNS TTETSSTLEV RTQGPILTAS GKNPVMELNE KRRGLKYELI SETGGSHDKR FVMEVEVDGQ KFRGAGPNKK VAKASAALAA LEKLFSGPNA ANNKKKKIIP QAKGVVNTAV SAAVQAVRGR GRGTLTRGAF VGATAAPGYI

image



Background

Alternative Names Spermatid perinuclear RNA-binding protein

Note

Product Name: Recombinant Human STRBP Protein
Catalog #: PEH9005



For Research Use Only , Not for Diagnostic Use