


Product Name: Recombinant Human SPSB2 Protein
Catalog #: PEH9007




Summary

Name	SPSB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99619
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSTPTPQALY PDLSCPEGLE ELLSAPPPDL GAQRRHGWNP KDCSENIEVK EGGLYFERRP VAQSTDGARG KRGYSRGLHA WEISWPLEQR GTHAVVGVAT ALAPLQTDHY AALLGSNSES WGWDIGRGKL YHQSKGPGAP QYPAGTQGEQ LEVPERLLVV LDMEEGTLGY AIGGTYLGPA FRGLKGRTLY PAVSAVWGQC QVRIRYLGGER

image



SDS



WB

Background

Alternative Names SPRY domain-containing SOCS box protein 2 (SSB-2) (Gene-rich cluster protein C9)

Note

Product Name: Recombinant Human SPSB2 Protein
Catalog #: PEH9007



For Research Use Only , Not for Diagnostic Use