


Product Name: Recombinant Human SPSB4 Protein
Catalog #: PEH9008



Summary

Name	SPSB4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96A44
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAKRELRGAE PGRPARLDQL LDMPAAGLAV QLRHAWNPED RSLNVFVKDD DRLTFHRHPV AQSTDGIRGK VGHARGLHAW QINWPARQRG THAVVGVATA RAPLHSVGYT ALVGSDAESW GVDLGRSRLY HDGKNQPGVA YPAFLGPDEA FALPDSLLVV LDMDEGTLSE IVDGQYLGVA FRGLKGGKLY PVVSAVWGHC EVTMRYINGL

image



Background

Alternative Names SPRY domain-containing SOCS box protein 4 (SSB-4)

Note

Product Name: Recombinant Human SPSB4 Protein
Catalog #: PEH9008



For Research Use Only , Not for Diagnostic Use