


Product Name: Recombinant Human SRRD Protein
Catalog #: PEH9009



Summary

Name	SRR1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UH36
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEAEFESDSG VVLRRIWEAE KDLFISDFWS SALETINRCL TKHLEQLKAP VGTLSDIFGN LHLDSLPEES DVATDSIPRE ILVTGTCHLK CVCYGIGNFA TCIVARNQLT FLLLLLEKCQ IPRSHCWVYD PLFSQLEIEV LNTLGVTVLS ENEEGKRSIR GEPTIFYMLH CGTALYNNLL WSNWSVDALS KMVIIGNSFK GLEERLLARI LQKNYPYIAK ILKGLEEF PQT SQYMDIF NDTSVHWFPV QKLEQLSIDI WEFREEPDYQ

image



Background

Alternative Names SRR1-like protein (SRR1 domain-containing protein)

Product Name: Recombinant Human SRRD Protein
Catalog #: PEH9009



Note

For Research Use Only , Not for Diagnostic Use