

Product Name: Recombinant Human STARD5 Protein
Catalog #: PEH9010



Summary

Name	STAR5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9NSY2
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDPALAAQMS EAVA EKMLQY RRDTAGWKIC REGNGVSVSW RPSVEFPGNL YRGE GIVYGT LEEVWDCVKP AVGGLRVKWD ENVTGF EIIQ SITDTLCVSR TSTPSAAMKL ISPRDFVDLV LVKRYEDGTI SSNATHVEHP LCPPKPGFVR GFNHPCGCFC EPLPGEPTKT NLVTFFHTDL SGYLPQNVVD SFFPRSMTRF YANLQKAVKQ FHE

image



Background

Alternative Names	StAR-related lipid transfer protein 5 (START domain-containing protein 5) (StARD5)
--------------------------	--

Note

Product Name: Recombinant Human STARD5 Protein
Catalog #: PEH9010



For Research Use Only , Not for Diagnostic Use