


Product Name: Recombinant Human SULT1C3 Protein
Catalog #: PEH9012



Summary

Name	ST1C3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6IMI6
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FYEDIKKDPK REIEKILKFL EKDISEEILN KIIYHTSFDV MKQNPMNTNYT TLPTSIMDHS ISPFMRKGMP GDWKNYFTVA QNEEFDKDYQ

image



Background

Alternative Names Sulfotransferase 1C3 (ST1C3) (EC 2.8.2.1) (EC 2.8.2.2)

Note

For Research Use Only , Not for Diagnostic Use