


Product Name: Recombinant Human SULT6B1 Protein
Catalog #: PEH9013



Summary

Name	ST6B1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6IMI4
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EFPVLECGDS EKYQRMKGFP SPRILATHLH YDKLPGSIFE NKAKILVIFR NPKDTAVSFL HFNHNDVPDIP SYGSWDEFFR QFMKGQVSWG RYFDFAINWN KHLGDGDNVVF ILYEDLKENL AAGIKQIAEF LGFFLTGEQI QTISVQSTFQ AMRAKSQDTH GAVGPFLFRK GEVGDWKNLF SEIQNQEMDE

image



Background

Alternative Names Sulfotransferase 6B1 (ST6B1) (Thyroxine sulfotransferase) (EC 2.8.2.n2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SULT6B1 Protein
Catalog #: PEH9013

