

Product Name: Recombinant Human SYNM Protein
Catalog #: PEH9015



Summary

Name	SYNEM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15061
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QHRRTKQPQE NTTHVEEVTE AGDSEGEQSY FVSTPDEHPG GHDRDDGSVY GQIHIEEEST IRYSWQDEIV QGTRRRRTQKD GAVGEKVVKP LDVPAPSLEG DLGSTHWKEQ ARSGEFHAEP TVIEKEIKIP HEFHTSMKGI SSKEPRQQLV EVIGQLEETL PERMREELSA LTREGQGGPG SVSVDVKKVQ GAGGSSVTLV

image



Background

Alternative Names Synemin (Desmuslin)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SYNM Protein
Catalog #: PEH9015

