


Product Name: Recombinant Human TBKBP1 Protein
Catalog #: PEH9016



Summary

Name	TBKB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A7MCY6
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ILTQEALGPS EVWLDSPGDP SLGGDMCSAS HFALITAYGD IKERLGGLER ENATLRRRLK VYEIKYPLIS DFGEEHGFSL YEIKDGSLL EKVSLQQRL NQFQHELQKN KEQEEQLGEM IQAYEKLCEV KSDLETRE MRALVETHLR QICGLEQQLR

image



Background

Alternative Names TANK-binding kinase 1-binding protein 1 (TBK1-binding protein 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TBKBP1 Protein
Catalog #: PEH9016

