


Product Name: Recombinant Human TTBK1 Protein
Catalog #: PEH9017



Summary

Name	TTBK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TCY1
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSTSTPPQQN TRQTAAMFGV VNVTPVPGDL LRENTEDVLQ GEHLSDQENA PPILPGRPSE GLGSPHLVP HPGGPEAEVW EETDVNRNKL RINIGKSPCV EEEQSRGMGV PSSPVRAPPD SPTTPVRSLR YRRVNSPESE RLSTADGRVE LPERRSRMDL PGSPSRQACS SQPAQMLSVD TGHADRQASG RMDVSASVEQ EALSNAFRSV PLAEEDFDS KEWVIIDKET

image



Background

Alternative Names Tau-tubulin kinase 1 (EC 2.7.11.1) (Brain-derived tau kinase)

Note

Product Name: Recombinant Human TTBK1 Protein
Catalog #: PEH9017



For Research Use Only , Not for Diagnostic Use