


Product Name: Recombinant Human TBC1D10A Protein
Catalog #: PEH9019



Summary

Name	TB10A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXI6
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AQKEQRKQMK GRGQLEKPPA PNQAMVVAAA GDACPPQHVP PKDSAPKDSA PQDLAPQVSA HHRSQESLTS

image



Background

Alternative Names	TBC1 domain family member 10A (EBP50-PDX interactor of 64 kDa) (EPI64 protein) (Rab27A-GAP-alpha)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use