


Product Name: Recombinant Human TEK2 Protein
Catalog #: PEH9020



Summary

Name	TEK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UIF3
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MATLSVKPSR RFQLPDWHTN SYLLSTNAQL QRDASHQIRQ EARVLRNETN NQTIWDEHDN RTRLVERIDT VNRWKEMLDK CLTDLDAEID ALTQMKESAE QNLQAKNLPL DVAIECLTLR ESRRDIDVVK DPVEDELHKE VEVIEATKKA LQQKVSQAFE QLCLLQEVQQ QLNSDHRGKM ETLEIDRGCL SLNLRSPNIS LKVDPTRVPD GSTTLQQWDD FSRFNKDRAE AEMKAATELR EATALTIAET NNELEAQRVA TEFAFRKRLR EMEKVYSELK WQEKNTLEEI AELQEDIRHL EEDLRTKLLS LKLSHTRLEA RTYRPNVELC RDQAQYGLTD EVHQLEATIA ALKQKLAQAQ DALDALCKHL ARLOADIACK ANSMLLDTKC MDTRRKLTVP AERFVPEVDT FTRTTNSTLS PLKSCQLELA

image



Background

Product Name: Recombinant Human TEKT2 Protein
Catalog #: PEH9020



Alternative Names

Tektin-2 (Tektin-t) (Testicular tektin) (Testicular tektin B1-like protein)
(TEKTB1) (Tektin-B1)

Note

For Research Use Only , Not for Diagnostic Use