


**Product Name: Recombinant Human TEKT2 Protein**  
**Catalog #: PEH9021**



## Summary

<b>Name</b>	TEKT2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UIF3
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	43 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SYLLSTNAQL QRDASHQIRQ EARVLRNETN NQTIWDEHDN RTRLVERIDT VNRWKEMLDK CLTDLDAEID ALTQMKESAE QNLQAKNLPL DVAIECLTLR ESRRDIDVVK DPVEDELHKE VEVIEATKKA LQQKVSQAFE QLCLLQEVQQ QLNSDHRGKM ETLEIDRGCL SLNLRSPNIS LKVDPTRVPD GSTTLQQWDD FSRFNKDRAE AEMKAATELR EATALTIAET NNELEAQRVA TEFAFRKRLR EMEKVYSELK WQEKNTLEEI AELQEDIRHL EEDLRTKLLS LKLSHTRLEA RTYRPNVELC RDQAQYGLTD EVHQLEATIA ALKQKLAQAA DALDALCKHL

## image



SDS

## Background

**Alternative Names** Tektin-2 (Tektin-t) (Testicular tektin) (Testicular tektin B1-like protein) (TEKTB1) (Tektin-B1)

**Product Name: Recombinant Human TEKT2 Protein**  
**Catalog #: PEH9021**

---



**Note**

For Research Use Only , Not for Diagnostic Use