



Product Name: Recombinant Human C9orf78 Protein
Catalog #: PEH9024



Summary

Name	TLS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZ63
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETTLVDDPFQ MKTGGMVDMK KLKERGKDKI SEEDLHLGT SFSAETNRRD EDADMMKYIE TELKKRKGIV EHEEQKVKPK NAEDCLYELP ENIRVSSAKK TEEMLSNQML SGIPEVDLGI DAKIKNIIST EDAKARLLAE QQNKKKDSET SFVPTNMAVN YVQHNRFYHE ELNAPIRRNK

image



Background

Alternative Names	Telomere length and silencing protein 1 homolog (Hepatocellular carcinoma-associated antigen 59)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C9orf78 Protein
Catalog #: PEH9024

