


Product Name: Recombinant Human TENT5A Protein
Catalog #: PEH9025



Summary

Name	TET5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96IP4
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEDELACSPY IPLGGDFGGG DFGGGDFGGG DFGGGGSFGG HCLDYCESPT AHCNVLNWEQ VQRLDGILSE TIPIHGRGNF PTLELQPSLI VKVVRRLAE KRIGVRDURL NGSAAASHVLH QDSGLGYKDL DLIFCADLRG EGEFQTVKDV VLDCLLDFLP EGVNKEKITP LTLKEAYVQK MVKVCNDSDR WSLISLSNNS GKNVELKFVD SLRRQFEFSV DSFQIKLDSL LLFYECSENP MTETFHPTII GESVYGDFQE AFDHLCNKII ATRNP EEIRG GLLKYCNLL VRGFRPASDE

image



Background

Alternative Names	Terminal nucleotidyltransferase 5A (EC 2.7.7.-) (HBV X-transactivated gene 11 protein) (HBV XAg-transactivated protein 11)
--------------------------	--

Product Name: Recombinant Human TENT5A Protein
Catalog #: PEH9025



Note

For Research Use Only , Not for Diagnostic Use