


Summary

Name	TTC19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6DKK2
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGAAAEDGAD EAEAEIIQLL KRAKLSIMKD EPEEAELILH DALRLAYQTD NKKAITTYTD LMANLAFIRG QLENAEQLFK ATMSYLLGGG MKQEDNAIIE ISLKLASIYA AQNRQEFAVA GYFCISTLE EKIEREKELA EDIMSVEEKA NTHLLGMCL DACARYLLFS KQPSQAQRMV EKALQISEEI

image



Background

Alternative Names Tetratricopeptide repeat protein 19, mitochondrial (TPR repeat protein 19)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TTC19 Protein
Catalog #: PEH9027

