


Product Name: Recombinant Human TRAPPC3L Protein
Catalog #: PEH9033



Summary

Name	TPC3L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T215
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSRPAHRRPE YHKINKDLFV LTYGALVAQL CKDYEKDEDV NQYLDKMGYG IGTRLVEDFL ARSCVGRCHS YSEIIDIIAQ VAFKMYLGIT PSVTCNNSSK NEFSLILEKN PLVEFVEELP AGRSSLCYCN LLCGIIRGAL EMVHLAADVT FLQDRLKGDS VTEIGITFLK KRDEKKYRGK K

image



Background

Alternative Names Trafficking protein particle complex subunit 3-like protein (TRAPPC3-like protein) (BET3-like protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRAPPC3L Protein
Catalog #: PEH9033

