


Product Name: Recombinant Human MTERF3 Protein
Catalog #: PEH9037



Summary

Name	MTEF3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96E29
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EELDELPPLS PMQPISEEEA IQIADPPLP PASFTLRDYV DHSETLQKLV LLGVDLSKIE KHPEAANLLL RLD FEKDIKQ MLLFLKDVGI EDNQLGAFLT KNHAIFSEDL ENLKTRVAYL HSKNFSKADV

image



Background

Alternative Names	Transcription termination factor 3, mitochondrial (Mitochondrial transcription termination factor 3) (mTERF3) (mTERF domain-containing protein 1, mitochondrial)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MTERF3 Protein
Catalog #: PEH9037

