


Product Name: Recombinant Human TBL2 Protein
Catalog #: PEH9038



Summary

Name	TBL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4P3
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRVFKMTKRE DGGYTFTATP EDFPKKHKAP VIDIGIANTG KFIMTASSDT TVLIWSLKGQ VLSTINTNQM NNTHAAVSPC GRFVASCFT PDVKVWEVCF GKKGEFQEVV RAFELKGHSA AVHSFAFSND SRRMASVSKD GTWKLWDTDV EYKKKQDPYL LKTGRFEEAA GAAPCRLALS PNAQVLALAS GSSIHLYNTR RGEKEECFER VHGECIANLS

image



Background

Alternative Names	Transducin beta-like protein 2 (WS beta-transducin repeats protein) (WS-betaTRP) (Williams-Beuren syndrome chromosomal region 13 protein)
--------------------------	---

Note

Product Name: Recombinant Human TBL2 Protein
Catalog #: PEH9038



For Research Use Only , Not for Diagnostic Use