



Product Name: Recombinant Human TACO1 Protein
Catalog #: PEH9040



Summary

Name	TACO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSH4
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKDTYLLYEG RGPGGSSLLI EALSNSSHKC QADIRHILNK NGGVMVAVGAR HSFDKKGIV VEVEDREKKA VNLERALEMA IEAGAEDVKE TEDEEERNVF KFICDASSLH QVRKKLDSLGLCSVSCALEF

image



Background

Alternative Names	Translational activator of cytochrome c oxidase 1 (Coiled-coil domain-containing protein 44) (Translational activator of mitochondrially-encoded cytochrome c oxidase I)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TACO1 Protein
Catalog #: PEH9040

