


Product Name: Recombinant Human TRIM34 Protein
Catalog #: PEH9042



Summary

Name	TRI34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BYJ4
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEKLLLFCKE DRKVICWLCE RSQEHARGHHT VLTEEVFKEC QEKLQAVLKR LKKEEEAEK LEADIREECT SWKYQVQTER QRIQTEFDQL RSILNNEEQR ELQRLEEEEEK KTLDFKFAEAE

image



Background

Alternative Names	Tripartite motif-containing protein 34 (Interferon-responsive finger protein 1) (RING finger protein 21)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use