


Product Name: Recombinant Human TRIM40 Protein
Catalog #: PEH9043



Summary

Name	TRI40_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6P9F5
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKLQALQFQV DHGNHRLEAG PESQHQTREQ LGALPQQWLG QLEHMPAEAA RILDISRAVT QLRSLVIDLE RTAKELDTNT LKNAGDLLNR SAPQKLEVIY PQLEKGVSEL

image



Background

Alternative Names Tripartite motif-containing protein 40 (Probable E3 NEDD8-protein ligase) (RING finger protein 35)

Note

For Research Use Only , Not for Diagnostic Use