


Product Name: Recombinant Human TRIM67 Protein
Catalog #: PEH9046



Summary

Name	TRIM67_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q6ZTA4
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGLRGFQRNR LLEAIVQRYQ QGRGAVPGTS AAAAVAICQL CDRTTPEPAA TLCEQCDVLY CSACQLKCHP SRGPFKHLR VQPPPPPPP AEAASGPTGT AQGAPSGGGG CKSPGGAGAG ATGGSTARKF PTCPEHEMEN YSMYCVSCRT PVCYLCLEEG RHAKHEVKPL GAMWKQHKAQ LSQALNGVSD KAKEAKEFLV QLKNILQQIQ ENGLDYEACL VAQCDALVDA LTRQKAKLLT KVTKEREHL

image



Background

Alternative Names Tripartite motif-containing protein 67 (TRIM9-like protein)

Note

Product Name: Recombinant Human TRIM67 Protein
Catalog #: PEH9046



For Research Use Only , Not for Diagnostic Use