


Product Name: Recombinant Human TRMT10B Protein
Catalog #: PEH9049



Summary

Name	TM10B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PF06
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDWKLEGSTQ KVESPVLQGQ EGILEETGED GLPEGFQLLQ IDAEGECQEG EILATGSTAW CSKNVQRKQR HW EKIVAAKK SKRKQEKERR KANRAENPGI CPQH SKRFLR ALTKDKLLEA KHSGPRLCID LSMTHYMSKK ELSRLAGQIR RLYG SNKKAD RPFWICLTGF TTDSPLYEEC VRMNDGFSSY LLDITEEDCF SLFPLETLVY LTPDSEHALE DVDLNKVYIL GGLVDESIQK KVTFQKAREY SVKTARLPIQ EYMRNQNGK NYHSEILAIN QVFDILSTYL ETHNWPEALK KGVSSGKGYI LRNSVE

image



Background

Alternative Names	tRNA methyltransferase 10 homolog B (EC 2.1.1.221) (RNA (guanine-9)-methyltransferase domain-containing protein 3) (tRNA (guanine(9)-N(1))-
--------------------------	---

Product Name: Recombinant Human TRMT10B Protein
Catalog #: PEH9049



methytransferase TRMT10B)

Note

For Research Use Only , Not for Diagnostic Use