

Product Name: Recombinant Human C11orf80 Protein
Catalog #: PEH9052



Summary

Name	TO6BL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N6T0
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITQHQLHHC SCEVKQLTL EKKDSAQGTE DAPDNSSLEL LADTSGQAEN KRLKRGSPRI EEMRALRSAR APSPSEAAPR RPEATAAPLT

image



Background

Alternative Names Type 2 DNA topoisomerase 6 subunit B-like (Type 2 DNA topoisomerase VI subunit B-like) (TOPOVIBL)

Note

For Research Use Only , Not for Diagnostic Use