


Product Name: Recombinant Human USP20 Protein
Catalog #: PEH9054



Summary

Name	UBP20_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2K6
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKRRRIEIDT FIKLNKAFQA EESPGVIYCI SMQWFREWEA FVKGKDNEPP GPIDNSRIAQ VKGSGHVQLK QGADYGQISE ETWTYLNSLY GGGPEIAIRQ SVAQPLGPEN

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 20 (EC 3.4.19.12) (Deubiquitinating enzyme 20) (Ubiquitin thioesterase 20) (Ubiquitin-specific-processing protease 20) (VHL-interacting deubiquitinating enzyme 2) (hVDU2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USP20 Protein
Catalog #: PEH9054

