

Product Name: Recombinant Human UEVLD Protein
Catalog #: PEH9057



Summary

Name	UEVLD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IX04
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLGKYKFRDL TVEELRNVNV FFPFKYSMD TYVFKDSSQK DLLNFTGTIP VMYQGNTYNI PIRFWILDSH PFAPPICFLK PTANMGILVG KHVDAQGRIY LPYLQNWSPH KSVIVGLIKE MIAKFQEELP MYSLSSSDEA RQVDLLAYIA KITEGVSDTN SKSWANHENK TVNKITVVG GELGIAC TLA ISAKGIADRL VLLDLSEGTK GATMDLEIFN LPNVEISKDL SASAHSKVVI FTVNSLGSSQ SYLDVVQSNV DMFRALVPAL GHYSQHSVLL VASQPVEIMT YVTWKLSTFP ANRVIGIGCN LDSQRLQYII TNVLKAQTSG KEVWVIGEQQ EDKVLTWSGQ EEVVSHTSQV QLSNRAMELL RVKGQRWSV GLSVADMVDS IVNNKKKVHS VSALAKGYD INSEVFLSLP CILGTNGVSE VIKTTLKEDT VTEKLQSSAS

image



Background

Product Name: Recombinant Human UEVLD Protein
Catalog #: PEH9057



Alternative Names

Ubiquitin-conjugating enzyme E2 variant 3 (UEV-3) (EV and lactate/malate dehydrogenase domain-containing protein)

Note

For Research Use Only , Not for Diagnostic Use