


Product Name: Recombinant Human UCK1 Protein
Catalog #: PEH9062



Summary

Name	UCK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HA47
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRPFLIGVSG GTASGKSTVC EKIMELLGQN EVEQRQRKVV ILSQDRFYKV LTAEQKAKAL KGQYNFDHPD AFDNDLMHRT LKNIVEGKTV EVPTYDFVTH SRLPETTVVY PADVVLFEGL LVFYSQEIRD MFHLRLFVDT DSDVRLSRRV LRDVRGRDL EQILTQYTF VKPAFEFCL PTKKYADVII PRGVDNMVAI NLIVQHIQDI LNGDICKWHR GGSNGRSYKR TFSEPGDHPG

image



Background

Alternative Names Uridine-cytidine kinase 1 (UCK 1) (EC 2.7.1.48) (Cytidine monophosphokinase 1) (Uridine monophosphokinase 1)

Note

Product Name: Recombinant Human UCK1 Protein
Catalog #: PEH9062



For Research Use Only , Not for Diagnostic Use