


Product Name: Recombinant Human USP6NL Protein
Catalog #: PEH9063



Summary

Name	US6NL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92738
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RHEHSPHPQS RTGTPERAQP PRRKSVEEES KKLKDEADFQ RKLPSGPQDS SRQYNHAAAN QNSNATSNIR KEFVPKWNKP SDVSATERTA KYTMEGKGRA AHPALAVTVP GPAEVRVSNV RPKMKALDAE DGKRGSTASQ YDNVPGPELD SGASVEEALE RAYSQSPRHA LYPPSPRKHA EPSSSPSKVS

image



Background

Alternative Names USP6 N-terminal-like protein (Related to the N-terminus of tre) (RN-tre)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USP6NL Protein
Catalog #: PEH9063

