


Product Name: Recombinant Human MON1B Protein
Catalog #: PEH9064



Summary

Name	MON1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L1V2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKDKDQPPSP SPPQSEALS STSRLWSPAA PENSPTCSPE SSSGGQGGDP SDEEWSQRK HVFVLSEAGK PIYSRYGSVE ALSATMGVMT

image



Background

Alternative Names Vacuolar fusion protein MON1 homolog B (HSV-1 stimulation-related gene 1 protein) (HSV-1 stimulating-related protein)

Note

For Research Use Only , Not for Diagnostic Use