


Product Name: Recombinant Human ZC3H11A Protein
Catalog #: PEH9068



Summary

Name	ZC11A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75152
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKTLEEIKLE KALRVQQSSE SSTSSPSQHE ATPGARRLLR ITKRTGMKEE KNLQEGNEVD SQSSIRTEAK EASGETTGVD ITKIQVKRCE TMREKHMQKQ QEREKSVLTP LRGDVASCNT QVAEKPVLTA VPGITRHLTK RLPTKSSQKV EVETSGIGDS LLNVKCAAQT LEKRGKAKPK VNVKPSVVKV VSSPKLAPKR

image



Background

Alternative Names Zinc finger CCCH domain-containing protein 11A

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZC3H11A Protein
Catalog #: PEH9068

