


Product Name: Recombinant Human CCDC66 Protein
Catalog #: PEH9072



Summary

Name	CCD66_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A2RUB6
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVMANSKKTN FLRSMTALLD PAQIEERDRR RQKQLEHQKA ITAQVEEKRR KKQLEEEQRK KEEQEEELRL AQEREEMQKQ YEEDILKQKQ KEEIMTLKTN ELFQTMQRAQ ELAQRLKQEQ RIRELAQKGH DTSRLIKNLG VDTIQMEYNA SNISNSRHDS DEISGKMNTY MNSTTSKKDT GVQTDDLNIG IFTNAESHCG SLMERDITNC SSPEISAEI GQFSTKKNKQ

image



Background

Alternative Names Coiled-coil domain-containing protein 66

Note

Product Name: Recombinant Human CCDC66 Protein
Catalog #: PEH9072



For Research Use Only , Not for Diagnostic Use