


**Product Name: Recombinant Human NME4 Protein**  
**Catalog #: PEH9080**



## Summary

<b>Name</b>	NDKM_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O00746
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	18 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GGPSWTRERT LVAVKPDGVQ RRLVGDVIQR FERRGFTLVG MKMLQAPESV LAEHYQDLRR KPFYPALIRY MSSGPVVAMV WEGYNVVRAS RAMIGHTDSA EAAPGTIRGD FSVHISRNV I HASDSVEGAQ REIQLWFAQS ELVSWADGGQ

## image



## Background

<b>Alternative Names</b>	Nucleoside diphosphate kinase, mitochondrial (NDK) (NDP kinase, mitochondrial) (EC 2.7.4.6) (Nucleoside diphosphate kinase D) (NDPKD) (nm23-H4)
--------------------------	---

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human NME4 Protein**  
**Catalog #: PEH9080**

---

