


Product Name: Recombinant Human CDK2AP1 Protein
Catalog #: PEH9081



Summary

Name	CDKA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14519
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSSQYRQLLS DYGPPSLGYT QGTGNSQVPQ SKYAELLAI EELGKEIRPT YAGSKSAMER LKRGIIHARG LVRECLAETE

image



Background

Alternative Names Cyclin-dependent kinase 2-associated protein 1 (CDK2-associated protein 1) (Deleted in oral cancer 1) (DOC-1) (Putative oral cancer suppressor)

Note

For Research Use Only , Not for Diagnostic Use