


Product Name: Recombinant Human FCHO1 Protein
Catalog #: PEH9082



Summary

Name	FCHO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14526
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKNHGFVLY HSVKQGPIST KELADFIRER ATIEETYSKA MAKLSKLASN GTPMGTFAPL WEVFRVSSDK LALCHLELTR KLQDLIKDVL RYGEEQLKTH KKCKEEVVST LDAVQVLSGV SLLPKSREN YLNRCMDQER LRRESTSQKE MDKAETKTKK AAESLRRSVE KYNSARADFE QKMLDSALRF QAMEETHLRH MKALLGSYAH SVEDTHVQIG QVHEEFKQNI ENVSVEMLLR KFAESKGTGR EKPGPLDFA YSAAALQEAM KRLRGAKAFR LPGLSRRERE PEPPAAVDFL EPDSGTCPEV

image



Background

Alternative Names F-BAR domain only protein 1

Product Name: Recombinant Human FCHO1 Protein
Catalog #: PEH9082



Note

For Research Use Only , Not for Diagnostic Use