


Product Name: Recombinant Human GSTA4 Protein
Catalog #: PEH9086



Summary

Name	GSTA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15217
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NGRGRMESVR WVLAAAGVEF DEEFLETKEQ LYKLQDGNHL LFQQVPMVEI DGMKLVQTRS ILHYIADKHN LFGKLNKERT LIDMYVEGTL DLLELLIMHP FLKPDDQQKE VVNMAQKAI RYFPVFEKIL RGHGQSFLVG NQLSLADVIL LQTILALEEK IPNILSAFPF LQEYTVKLSN IPTIKRFLEP GSKKKPPPDE

image



Background

Alternative Names Glutathione S-transferase A4 (EC 2.5.1.18) (GST class-alpha member 4)
(Glutathione S-transferase A4-4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GSTA4 Protein
Catalog #: PEH9086

