


Product Name: Recombinant Human CASC3 Protein
Catalog #: PEH9087



Summary

Name	CASC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15234
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DAVNSSTKEE KGEEKPDTKS TVTGERQSGD GQESTEPVEN KVGKKGPKHL DDEDEDRKNPA YIPRKGLFFE HDLRGQTQEE EVRPKGRQRK LWKDEGRWEH DKFREDEQAP KSRQELIALY

image



Background

Alternative Names Protein CASC3 (Cancer susceptibility candidate gene 3 protein) (Metastatic lymph node gene 51 protein) (MLN 51) (Protein barentsz) (Btz)

Note

For Research Use Only , Not for Diagnostic Use