


Product Name: Recombinant Human MAGEB2 Protein
Catalog #: PEH9089



Summary

Name	MAGB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15479
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AREKRRKARD ETRGLNVPQV TEAEEEEAPC CSSSVSGGAA SSSPAAGIPQ EPQRAPTTAA AAAAGVSSTK SKKGAKSHQG EKNASSSQAS

image



Background

Alternative Names	Melanoma-associated antigen B2 (Cancer/testis antigen 3.2) (CT3.2) (DSS-AHC critical interval MAGE superfamily 6) (DAM6) (MAGE XP-2 antigen) (MAGE-B2 antigen)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use