


Product Name: Recombinant Human HSPA12A Protein
Catalog #: PEH9090



Summary

Name	HS12A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43301
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPRETAP TSA YSSPARSLGD TGITPLSPSH IVNDTDSNVS EQQSFLVVVA VDFGTTSSGY AYSFTKEPEC IHVMRRWEGG DPGVSNQKTP TTILLTPERK FHSFGYAARD FYHDLDPNEA KQWLYLEKFK MKLHTTGDLT MDTDLTAANG KKVKALEIFA YALQYFKEQA LKELSDQAGS EFENS DVRWV ITVPAIWKQP AKQFMRQAAY QAGLSPENS EQLIALEPE AASIYCRKLR LHQMIELSSK AAVNGYSGSD TVGAGFTQAK EHRRNRQSR TFLVENVIGE IWSELEEGDK YVVVDSGGGT VDLTVHQIRL PEGHLKELYK ATGGPYGSLG VDYEFKLLY KIFGEDFIEQ FKIKRPAAWV

image



Background

Alternative Names	Heat shock 70 kDa protein 12A
--------------------------	-------------------------------

Product Name: Recombinant Human HSPA12A Protein
Catalog #: PEH9090



Note

For Research Use Only , Not for Diagnostic Use