


**Product Name: Recombinant Human XPOT Protein**  
**Catalog #: PEH9092**



## Summary

<b>Name</b>	XPOT_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O43592
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	48 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SPDAWQVCAE ALAQRYSDD HVKFFCFQVL EHQVKYKYSE LTTVQQQLIR ETLISWLQAQ MLNPQPEKTF IRNKAQVFA LLFVTEYLTK WPKFFFDILS VVDLNPRGVD LYLRILMAID SELVDRDVVH TSEEARRNTL IKDTMREQCI PNLVESWYQI LQNYQFTNSE VTCQCLEVVG AYVSWIDLSL IANDRFINML LGHMSIEVLR EEACDCLFEV VNKGMDPVDK MKLVESLCQV LQSAGFFSID QEEDVDFLAR FSKLVNGMGQ SLIVSWSKLI KNGDIKNAQE ALQAIETKVA LMLQLLIHED DDISSNIIGF CYDYHLILKQ LTVLSDQQKA NVEAIMLAVM KKLTYDEEYN FENEGEDEAM FVEYRKQLKL LLDRLAQVSP

## image



## Background

**Alternative Names** Exportin-T (Exportin(tRNA)) (tRNA exportin)

**Product Name: Recombinant Human XPOT Protein**  
**Catalog #: PEH9092**

---



**Note**

For Research Use Only , Not for Diagnostic Use