


Product Name: Recombinant Human ZNRD2 Protein
Catalog #: PEH9095



Summary

Name	ZNRD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60232
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FSWEPTEAE TKVLQARRER QDRISRLMGD YLLRGYRMLG ETCADCGTIL LQDKQRKIYC VACQELDSDV DKDNPALNAQ AALSQAREHQ LASASEPLG SRPAPQPPVP RPEHCEGAAA GLKAAQGPPA PAVPPNTDVM ACTQTALLQK LTWASAEELGS STSLETSIQL

image



Background

Alternative Names	Protein ZNRD2 (Autoantigen p27) (Sjogren syndrome/scleroderma autoantigen 1) (Zinc ribbon domain-containing protein 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNRD2 Protein
Catalog #: PEH9095

