


Product Name: Recombinant Human SRGAP2 Protein
Catalog #: PEH9097



Summary

Name	SRGP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75044
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IEVISEPPEE KVTARAGASC PSGGHVADIY LANINKQRKR PESGSIRKTF RSDSHGLSSS LTDSSSPGVG ASCRPSSQPI MSQSLPKEGP DKCSISGHGS LNSISRHSSL KNRLDSPQIR KTATAGRSKS FNNHRPMDPE VIAQDIEATM NSALNELREL

image



Background

Alternative Names	SLIT-ROBO Rho GTPase-activating protein 2 (srGAP2) (Formin-binding protein 2) (Rho GTPase-activating protein 34)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SRGAP2 Protein
Catalog #: PEH9097

