


Product Name: Recombinant Human PEX14 Protein
Catalog #: PEH9100



Summary

Name	PEX14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75381
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASVQELLIQQ QKIQELAHE LAAAKATTST NWILESQNIN ELKSEINSLK GLLLNRRQFP PSPSAPKIPS WQIPVKSPSP SSPAAVNHHS SSDISPVSNE STSSSPGKEG HSPEGSTVTY HLLGPQEEGE GVVDVKGQVR MEVQGEEER EDKEDEEED

image



Background

Alternative Names Peroxisomal membrane protein PEX14 (PTS1 receptor-docking protein) (Peroxin-14) (Peroxisomal membrane anchor protein PEX14)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PEX14 Protein
Catalog #: PEH9100

