

Product Name: Recombinant Human AKAP3 Protein
Catalog #: PEH9103



Summary

Name	AKAP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75969
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEHVRKAQDK AESYLSIMK GMGDPKNNRV NFAMKSETKL REKMYSEPKS EEETCAKTLG EHIKEGLTL WHKTQQKECK

image



Background

Alternative Names	A-kinase anchor protein 3 (AKAP-3) (A-kinase anchor protein 110 kDa) (AKAP 110) (Cancer/testis antigen 82) (CT82) (Fibrous sheath protein of 95 kDa) (FSP95) (Fibrousheathin I) (Fibrousheathin-1) (Protein kinase A-anchoring protein 3) (PRKA3) (Sperm oocyte-binding protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AKAP3 Protein
Catalog #: PEH9103

