


Product Name: Recombinant Human BAIAP3 protein
Catalog #: PEH9104



Summary

Name	BAIP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94812
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQDPGSASAD PQEPATGAWK PGDGVEFFAH MRLMLKKGEG RQGLPCLEVP LRSGSPAPPE PVDPSLGLRA LAPEEVEMLY EEALYTVLYR AGTMGPDQVD DEEALLSYLQ QVFGTSLEEH TEAIERVRKA KAPTYALKVS VMRAKNLLAK DPNGFSDPYC MLGILPASDA TREPRAQKEQ RFGFRKGSKR GGPLPAKCIQ VTEVKSSTLN PVWKEHFLFE IEDVSTDQLH LDIWDHDDDV SLVEACRKLN EVIGLKGMMGR YFKQIVKSAR ANGTAIPTED HTDDFLGCLN

image



Background

Alternative Names	BAI1-associated protein 3 (BAP3) (Brain-specific angiogenesis inhibitor I-associated protein 3)
--------------------------	---

Product Name: Recombinant Human BAIAP3 protein
Catalog #: PEH9104



Note

For Research Use Only , Not for Diagnostic Use