


Product Name: Recombinant Human GAS8 Protein
Catalog #: PEH9108



Summary

Name	DRC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95995
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITRMRNDFER QVREIEAKYD KKMKMLRDEL DLRRKTELHE VEERKNGQIH TLMQRHEEAF TDIKNYNDI TLNNLALINS LKEQMEDMRK KEDHLEREMA EVSGQNKRLA DPLQKAREEM SEMQQLANY

image



Background

Alternative Names Dynein regulatory complex subunit 4 (Growth arrest-specific protein 11) (GAS-11) (Growth arrest-specific protein 8) (GAS-8)

Note

For Research Use Only , Not for Diagnostic Use